

Setting up your CPAP machine

Ph: 1800 799 950

E: info@sleepright.com.au

W: sleepright.com.au

The Equipment:

The 5 essential components of CPAP are

- CPAP device - Humidifier - Power Supply - Tubing - CPAP Mask

How to set up your CPAP machine:

1. Place the machine on a stable, flat surface at head height or below



2. Plug the power supply unit into the wall and into the back of the CPAP machine. Ensure that the cord is connected properly at the transformer



3. If the machine has a humidifier, fill it with water (make sure you do not overfill it by taking note of the maximum fill line)

4. Connect the tubing to the mask & machine. Heated tubes have a slot for the power chip to be connected



5. Press the on or power button to begin and stop therapy



***If you have any questions, please get in touch with us at info@sleepright.com.au or call us on 1800 799 950.